
pH Testing

pH paper is a key tool to test your body's overall acid/alkaline balance by testing the pH of your first morning urine (the first urine after 5 AM). This pH roll can test a range from 5.5 to 8.0 in steps of 0.2 pH units. One roll will last for nearly two hundred tests. A convenient dispenser allows you to tear off a small strip of pH paper (1 inch or less) to dip into urine to determine its pH.



pH Paper With pH Color Chart (0.4 oz)

This convenient and simple tool is perfect for monitoring your progress in helping to restore your body's acid/alkaline balance. You can easily check either urine or saliva pH.



Health Benefits

pH management is the cornerstone of health maintenance. The pH Paper roll offers key advantages:

- Easily track your progress.
- Conveniently fits in a shirt pocket or purse.
- Inexpensive
- One roll lasts several months

Synergistic Considerations

For maximum replenishment of pH reserve (if tested numbers fall in the yellow or blue zones), use the **pH Trio** (Coral Legend, Quantum Aloe and Quantum Cod Liver Oil) until the first morning urine pH falls consistently in the 6.4 to 7.0 range.

Directions

Withdraw and tear off approximately 1 inch or less of pH paper and immerse in a cup of freshly voided first morning urine. Remove the strip, shake off any excess urine and then promptly read the color; match the specimen's color with the colors on the pH dispenser chart to determine the corresponding pH.

Urine pH Tracking

Start Date_	
Name	

Please record your urine pH values each day on this chart.

AM pH reading is taken from testing the first morning urine after 5 AM. PM pH reading is taken from testing urine taken just before dinner meal (approx. 5-6 PM).

PIV																
AM																
	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
	Day															
PM																
AM																
	1	2	3	4	5	9	7	8	6	10	11	12	13	14	15	
	Day															

Urine pH Tracking

Date_		
NameStartStartStartStart Please record your urine pH values each day on this o	Vame Start Date	Please record your urine pH values each day on this chart.

AM pH reading is taken from testing the first morning urine after 5 AM. PM pH reading is taken from testing urine taken just before dinner meal (approx. 5-6 PM).

PM

AM																
	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
	Day															
PM																
AM																
	1	2	3	4	5	9	7	8	6	10	11	12	13	14	15	
	Day															